

StitchKit: multiple configurations for maximum flexibility

StitchKit

Pre-loaded with 6 sutures




Built-in Sharps Container

StitchKit comes pre-loaded in various suture combinations for maximum ease-of-use. Suture options include barbed, absorbable, non-absorbable, monofilament, and braided

StitchKit PARK

Load up to 6 sutures



Built-in Sharps Container

StitchKit PARK comes unloaded for maximum flexibility, allowing up to 6 sutures of any type to be loaded on the back table and deployed into the surgical field

Sharps Container



All StitchKit configurations have an on-board sharps container to ensure that needles are secure throughout the case



Stitch Kit
Videos



Stitch Kit
Configurations

Pre-loaded Configurations

	SKU	Quantity	Size	Type	Length (approx.)	Needle
Urogynecology	SK-002	6	2-0	PTFE	8"	V-20
		2	2-0	PTFE	8"	V-20
	SK-004	3	3-0	V-LOC 180	9"	V-20
		1	2-0	MAXON	9"	V-20
	SK-106	2	2-0	PTFE	8"	V-20
		4	2-0	MAXON	8"	V-20
General, Bariatric, Hernia	SK-103	1	2-0	POLYSORB	8"	V-20
		2	3-0	V-LOC 90	9"	V-20
		3	3-0	V-LOC 180	9"	V-20
	SK-104	3	2-0	SOFSILK	8"	V-20
		1	3-0	V-LOC 90	12"	V-20
		2	2-0	POLYSORB	8"	V-20
	SK-107	3	2-0	POLYSORB	8"	V-20
		3	3-0	V-LOC 180	9"	V-20
	SK-108	3	2-0	Polysorb	9"	V-20
		3	2-0	V-LOC 180	6"	GS-22
	SK-109	1	3-0	POLYSORB	8"	V-20
		3	3-0	V-LOC 90	6"	V-20
		2	3-0	V-LOC 180	6"	V-20
	SK-111	2	3-0	TICRON	8"	V-20
		1	3-0	POLYSORB	8"	V-20
		3	3-0	V-LOC 180	6"	V-20
	SK-211	6	2-0	TICRON	9"	V-20
Urology	SK-105	2	3-0	V-LOC 90	6"	V-20
		4	0	POLYSORB	8"	V-20
	SK-110	2	0	POLYSORB	8"	V-20
		2	2-0	POLYSORB	8"	CV-23
		1	2-0	BIOSYN (undyed)	8"	GU-46
		1	2-0	BIOSYN (dyed)	8"	GU-46
All Specialties	SK-101	Unloaded Configuration Ability to load any combination of suture in the Operating Room				

To order: 1.973.765.6256

Origamisurgical.com

Copyright 2025 Origami Surgical Inc, All Rights Reserved. StitchKit™ is a registered trademark of Origami Surgical Inc

Origami Surgical and the Origami Surgical logo are trademarks of Origami Surgical Inc.

This device is covered by US Patent 8,418,851 & 6,986,780.
 Other patents pending. V-loc™, Polysorb™, Sofsilk™, Biosyn™,
 Maxon™, and Ticon™ are trademarks of Medtronic

